

## Product Datasheet

# MAPKAPK2 (phospho-Thr334) antibody (HRP) orb560845

<b>Description</b>	Rabbit polyclonal antivosy against MAPKAPK2 (phospho-Thr334) conjugated to HRP
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	IHC-P
<b>Immunogen</b>	KLH conjugated synthetic peptide derived from the region around the phosphorylation site of T human MAPKAPK2 (phospho-Thr334)
<b>Target</b>	MAPKAPK2 (phospho-Thr334)
<b>Form/Appearance</b>	Liquid: 0.01M TBS (pH7.4), 0.5% BSA, 0.015% ProClin™ 300 and 25% glycerol
<b>Concentration</b>	- 100 µg (in 200 µl): 0.5 mg/ml
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Affinity purified by Protein A
<b>MW</b>	46 kDa
<b>Uniprot ID</b>	<b>P49137</b> (Human); <b>P49138</b> (Mouse)
<b>Entrez</b>	<b>9261</b>
<b>Dilution Range</b>	IHC-P: 1:50-200
<b>Expiration Date</b>	12 months from date of receipt.