

## Product Datasheet

### ADMR antibody (HRP) orb560870

<b>Description</b>	Rabbit polyclonal antivosy against ADMR conjugated to HRP
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Mouse, Rat
<b>Predicted Reactivity</b>	Mouse
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	ELISA, IHC-P
<b>Immunogen</b>	KLH conjugated synthetic peptide derived from mouse ADMR
<b>Target</b>	ADMR
<b>Form/Appearance</b>	Liquid: 0.01M TBS (pH7.4), 0.5% BSA, 0.015% ProClin™ 300 and 25% glycerol
<b>Concentration</b>	- 100 µg (in 200 µl): 0.5 mg/ml
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Affinity purified by Protein A
<b>MW</b>	45 kDa
<b>Uniprot ID</b>	<b>O15218</b> (Human); <b>P43142</b> (Mouse); <b>P31392</b> (Rat)
<b>Entrez</b>	<b>11318</b>
<b>Dilution Range</b>	ELISA: 1:100-1000, IHC-P: 1:50-200
<b>Expiration Date</b>	12 months from date of receipt.