

## Product Datasheet

### GABRD antibody (HRP) orb560874

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Rabbit polyclonal antivosy against GABRD conjugated to HRP   |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human, Mouse, Rat  |
| <b>Conjugation</b>         | HRP  |
| <b>Tested Applications</b> | ELISA, ICC, IHC-P  |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide derived from human GABRD  |
| <b>Target</b>              | GABRD  |
| <b>Form/Appearance</b>     | Liquid: 0.01M TBS (pH7.4), 0.5% BSA, 0.015% ProClin™ 300 and 25% glycerol                                      |
| <b>Concentration</b>       | - 100 µg (in 200 µl): 0.5 mg/ml  |
| <b>Storage</b>             | Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles. |
| <b>Note</b>                | For research use only.   |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Purity</b>              | Affinity purified by Protein A   |
| <b>MW</b>                  | 49 kDa   |
| <b>Uniprot ID</b>          | <b>O14764</b> (Human); <b>P22933</b> (Mouse); <b>P18506</b> (Rat)  |
| <b>Entrez</b>              | <b>2563</b>  |
| <b>Dilution Range</b>      | ELISA: 1:100-1000, IHC-P: 1:50-200, ICC: 1:50-200  |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |