

Product Datasheet

Cleaved PARP antibody (HRP) (orb560992)

Description	Mouse monoclonal antivosy against Cleaved PARP conjugated to HRP...
Species/Host	Mouse
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA, ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Cleaved PARP
Target	Cleaved PARP
Form/Appearance	Liquid: 0.01M TBS (pH7.4), 0.5% BSA, 0.015% ProClin™ 300 and 25% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Monoclonal
Purity	Affinity purified by Protein G
MW	89/116 kDa
Clone ID	1A1
Uniprot ID	P09874 (Human); P11103 (Mouse); P27008 (Rat)
Entrez	142
Dilution Range	WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200, ICC: 1:50-200
Expiration Date	12 months from date of receipt.