

Product Datasheet

STAT1 antibody (HRP) orb561020

| | |
|----------------------------|----------------------------------------------------------------------------------------------------------------|
| Description | Mouse monoclonal antibody against STAT1 conjugated to HRP |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | HRP |
| Tested Applications | WB |
| Immunogen | Recombinant human STAT1 |
| Target | STAT1 |
| Form/Appearance | Liquid: 0.01M TBS (pH7.4), 0.5% BSA, 0.015% ProClin™ 300 and 25% glycerol |
| Concentration | - 100 µg (in 200 µl): 0.5 mg/ml |
| Storage | Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles. |
| Note | For research use only. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | Affinity purified by Protein G |
| MW | 82 kDa |
| Clone ID | 5G8 |
| Uniprot ID | P42224 (Human); P42225 (Mouse) |
| Dilution Range | WB: 1:50-200 |
| Expiration Date | 12 months from date of receipt. |