

Product Datasheet

CELA1 antibody (HRP) orb561105

Description	Rabbit polyclonal antivosy against CELA1 conjugated to HRP
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	HRP
Tested Applications	ELISA, ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human CELA1
Target	CELA1
Form/Appearance	Liquid: 0.01M TBS (pH7.4), 0.5% BSA, 0.015% ProClin™ 300 and 25% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	28 kDa
Uniprot ID	Q9UNI1 (Human); P00773 (Rat)
Entrez	1990
Dilution Range	WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200, ICC: 1:50-200
Expiration Date	12 months from date of receipt.