

Product Datasheet

C1orf52 antibody (HRP) orb561232

Description	Rabbit polyclonal antivosy against C1orf52 conjugated to HRP
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	HRP
Tested Applications	ELISA, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human C1orf52
Target	C1orf52
Form/Appearance	Liquid: 0.01M TBS (pH7.4), 0.5% BSA, 0.015% ProClin™ 300 and 25% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	21 kDa
Uniprot ID	Q8N6N3 (Human)
Entrez	148423
Dilution Range	ELISA: 1:100-1000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.