

Product Datasheet

DDR1 antibody orb570379

Description	Rabbit polyclonal antibody to DDR1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FACS, IHC-P, WB
Immunogen	E.coli-derived human MCK10/NEP/DDR1 recombinant protein (Position: R341-A909).
Target	DDR1
Form/Appearance	Lyophilized: Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg NaN ₃ .
Storage	Store at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only.
Reconstitution	Add 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Isotype	IgG
Clonality	Polyclonal
Purity	Immunogen affinity purified.
Uniprot ID	Q08345
Entrez	780
Dilution Range	WB: 0.25-0.5 µg/ml, IHC-P: 0.5-1 µg/ml, FACS: 1-3 µg/1x10 ⁶ cells, D-ELISA: 0.1-0.5 µg/ml
Expiration Date	12 months from date of receipt.