

Product Datasheet

HTRA2 (Phospho-S142) antibody orb571001

Description	Rabbit polyclonal antibody to HTRA2 (Phospho-S142)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Custom conjugations	Unconjugated, Biotin, FITC, HRP, CF647, CF488A, CF594, CF405M
Tested Applications	ELISA, IHC-P
Immunogen	Synthesized peptide derived from human HtrA2 around the phosphorylation site of S142.
Target	HTRA2
Form/Appearance	Liquid: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Uniprot ID	O43464
Dilution Range	IHC-P: 1:100-1:300, ELISA: 1:40000
Expiration Date	12 months from date of receipt.