

## Product Datasheet

### NFKB2 (Phospho-S869) antibody orb571131

<b>Description</b>	Rabbit polyclonal antibody to NFKB2 (Phospho-S869)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Custom conjugations</b>	Unconjugated, Biotin, FITC, HRP, CF647, CF488A, CF594, CF405M
<b>Tested Applications</b>	ELISA, IHC-P, IP, WB
<b>Immunogen</b>	Synthesized peptide derived from human NFκB-p100 around the phosphorylation site of S869.
<b>Target</b>	NFKB2
<b>Form/Appearance</b>	Liquid: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Store at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Uniprot ID</b>	<b>Q00653</b>
<b>Dilution Range</b>	WB: 1:500-1:2000, IHC-P: 1:100-1:300, IP: 2-5 µg: mg lysate, ELISA: 1:10000
<b>Expiration Date</b>	12 months from date of receipt.