

Product Datasheet

STAT5A (Phospho-S726/731) antibody orb571316

Description	Rabbit polyclonal antibody to STAT5A (Phospho-S726/731)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Custom conjugations	Unconjugated, Biotin, FITC, HRP, CF647, CF488A, CF594, CF405M
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	Synthesized peptide derived from human Stat5a/b around the phosphorylation site of S726/731.
Target	STAT5A
Form/Appearance	Liquid: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Uniprot ID	P42229
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:100-1:300, IF: 1:200-1:1000, ELISA: 1:20000
Expiration Date	12 months from date of receipt.