

Product Datasheet

DNA Lambda Phage orb420124

Description	DNA Lambda Phage
Hazard Information	Non-Toxic
Target	DNA Lambda Phage
Form/Appearance	Liquid (sterile filtered)
Concentration	1450 µg/mL
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Expiration Date	12 months from date of receipt.
Application Notes	Purified Lambda phage DNA is in a solution in 10 mM Tris-Cl and 1 mM EDTA, pH 8.0. This product may be used for to perform activity assays for restriction enzymes, the preparation of DNA molecular weight standards, and other cloning procedures.