

Product Datasheet

Coating Stabilizer & Block Buf orb90373

Description	Coating Stabilizer & Blocking Buffer is a assay detection and miscellaneous reagent.
Hazard Information	Non-Toxic
Target	Coating Stabilizer & Blocking Buffer
Form/Appearance	pH 7.2
Storage	Store at -2 to -8°C.
Note	For research use only.
Expiration Date	12 months from date of receipt.
Application Notes	Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.