

Product Datasheet

Unmasking fluid T, TEC buffer (TRIS/ EDTA/Citrate solution) pH 8 orb90420

Description	Unmasking fluid T, TEC buffer (TRIS/ EDTA/Citrate solution) pH 8
Tested Applications	IHC-P
Hazard Information	Non-Toxic
Target	Unmasking fluid T, TEC buffer (TRIS/ EDTA/Citrate solution) pH 8
Form/Appearance	100 ml 10-times concentrated solution
Storage	Store at 4°C for short term and at -20°C for long term storage.
Note	For research use only.
Expiration Date	12 months from date of receipt.