

Product Datasheet

Mitochondrial Antigen orb179951

Description	Purified Autoimmune Reagents
Tested Applications	ELISA
Hazard Information	Non-Toxic
Target	Mitochondrial
Form/Appearance	Liquid in Buffer: 0.01M Tris containing 0.53M Sodium chloride, Preservative: 0.02% Sodium azide
Concentration	0.707 Units/ml
Storage	Store at -70°C. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only.
Purity	Affinity Purified
Expiration Date	12 months from date of receipt.