

Product Datasheet

2.5 mM dNTP Mix orb180488

Description	2.5 mM dNTP Mix consists of a solution of all four nucleotides (dATP, dCTP, dGTP, dTTP), each at a concentration of 2.5 mM. It is neutralized to pH 8.0 with NaOH, and supplied in purified water.
Hazard Information	Non-Toxic
Target	2.5 mM dNTP Mix
Form/Appearance	Supplied in purified water with NaOH and pH 8.0.
Concentration	2.5 mM
Storage	Store at -20°C
Note	For research use only.
Expiration Date	12 months from date of receipt.