

Product Datasheet

Beta galactosidase orb348524

Description	Beta Galactosidase
Hazard Information	Non-Toxic
Target	Beta galactosidase
Form/Appearance	Liquid: 0.05 M Potassium Phosphate, pH 7.8 10% (w/v) Sucrose
Concentration	1.0 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Purity	Beta Galactosidase was prepared from chromatographically pure beta-galactosidase (E.coli).
Expiration Date	12 months from date of receipt.
Application Notes	Beta Galactosidase is used as a control in beta-galactosidase based immunological assays.