

Product Datasheet

1.0 M Magnesium Chloride Depc Treated orb348571

| | |
|---------------------------|--|
| Description | 1.0 M Magnesium Chloride Depc Treated |
| Hazard Information | Non-Toxic |
| Target | 1.0 M Magnesium Chloride Depc Treated |
| Form/Appearance | Liquid (sterile filtered) |
| Concentration | 1.0 M |
| Storage | Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles. |
| Note | For research use only. |
| Purity | This product contains highly purified |
| Expiration Date | 12 months from date of receipt. |
| Application Notes | This product is a concentrated stock solution and should be diluted appropriately with distilled, deionized water or equivalent to its final working concentration. This buffer contains 1.0 M Magnesium Chloride and was meticulously prepared using ultra pure reagents dissolved in highly polished pharmaceutical grade deionized water treated with diethyl pyrocarbonate (DEPC). |