

Product Datasheet

Protract B Extraction Buffer orb348596

Description	Protract B Extraction Buffer
Tested Applications	WB
Hazard Information	Non-Toxic
Target	Protract B Extraction Buffer
Form/Appearance	Liquid (sterile filtered)
Concentration	1X
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Purity	This product contains highly purified
Expiration Date	12 months from date of receipt.
Application Notes	ProXtract B Buffer is suitable for protein extraction and easily scalable. Specific conditions should be optimized by the end user. See protocol.