

Product Datasheet

Borate Saline Buffered (Bsb) orb348617

Description	Borate Saline Buffered (Bsb)
Hazard Information	Non-Toxic
Target	Borate Saline Buffered (Bsb)
Form/Appearance	Liquid (sterile filtered)
Concentration	1X
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Purity	This product contains highly purified
Expiration Date	12 months from date of receipt.