

## Product Datasheet

### Bovine Serum Albumin 30% Solution orb348630

<b>Description</b>	Bovine Serum Albumin 30% Solution
<b>Tested Applications</b>	ELISA, ICC, IF, IHC-P, IP, WB
<b>Hazard Information</b>	Non-Toxic
<b>Target</b>	Bovine Serum Albumin 30% Solution
<b>Form/Appearance</b>	Liquid: 0.85% (w/v) Sodium Chloride 0.3% (v/v) Octanoic Acid
<b>Concentration</b>	30% (w/v)
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only.
<b>Purity</b>	This product contains highly purified
<b>Expiration Date</b>	12 months from date of receipt.