

Product Datasheet

Bovine Serum Albumin 30% Solution orb348632

Description	Bovine Serum Albumin 30% Solution
Tested Applications	ELISA, ICC, IF, IHC-P, IP, WB
Hazard Information	Non-Toxic
Target	Bovine Serum Albumin 30% Solution
Form/Appearance	Liquid: 0.85% (w/v) Sodium Chloride 0.3% (v/v) Octanoic Acid
Concentration	30% (w/v)
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Purity	This product contains highly purified
Expiration Date	12 months from date of receipt.