

## Product Datasheet

### Bovine Serum Albumin orb348633

<b>Description</b>	Bovine Serum Albumin
<b>Tested Applications</b>	ELISA, FACS, ICC, IF, IHC-P, IP, WB
<b>Hazard Information</b>	Non-Toxic
<b>Target</b>	Bovine Serum Albumin
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only.
<b>Purity</b>	Purity (%): 100% by Agarose Zone Electrophoresis Moisture (%):
<b>Expiration Date</b>	12 months from date of receipt.
<b>Application Notes</b>	This product is suitable for use in protease sensitive assays such as RIA, EIA and nucleic acid hybridization, use as a stabilizing agent for proteins and enzymes including dilute solutions of antibody, use as a blocking agent to reduce non-specific binding.