

## Product Datasheet

### HL-60 Whole Cell Lysate orb348702

<b>Description</b>	HL-60 Whole Cell Lysate
<b>Reactivity</b>	Human
<b>Tested Applications</b>	ChIP, IP, WB
<b>Hazard Information</b>	Non-Toxic
<b>Target</b>	HL-60 Whole Cell Lysate
<b>Form/Appearance</b>	Liquid: 1X RIPA Buffer with HALT Protease and Phosphatase Inhibitors
<b>Concentration</b>	4.0 mg/ml
<b>Storage</b>	Store vial at -70°C or COLDER. For extended storage, aliquot contents to minimize freeze/thaw cycles.
<b>Note</b>	For research use only.
<b>Expiration Date</b>	12 months from date of receipt.