

Product Datasheet

Streptavidin Atto 550 Conjugated orb348776

Description	Streptavidin Atto 550 Conjugated
Tested Applications	ICC, IF, WB
Hazard Information	Non-Toxic
Target	Streptavidin Atto 550
Form/Appearance	Lyophilized: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free
Concentration	1.1 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Purity	STREPTAVIDIN ATTO 550 Conjugated was prepared from chromatographically purified Streptavidin. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Streptavidin. No reaction was observed against anti-Avidin.
Expiration Date	12 months from date of receipt.
Application Notes	The emission spectra for this ATTO conjugate matches the principle output wavelengths of most common fluorescence instrumentation.