

Product Datasheet

1X Te, Ph 8.0 orb348559

Description	1X Te, Ph 8.0
Hazard Information	Non-Toxic
Target	1X Te, Ph 8.0
Form/Appearance	Liquid (sterile filtered)
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Purity	This product contains highly purified
Expiration Date	12 months from date of receipt.
Application Notes	This product is ready-to-use as a working 1X solution and requires no further dilution. This buffer consists of 10 mM Tris HCl and 1.0 mM EDTA at a pH of 8.0. Meticulously prepared using ultra pure reagents dissolved in DEPC treated pharmaceutical grade deionized water.