

Product Datasheet

5.0 M Potassium Acetate orb348561

Description	5.0 M Potassium Acetate
Hazard Information	Non-Toxic
Target	5.0 M Potassium Acetate
Form/Appearance	Liquid (sterile filtered)
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Purity	This product contains highly purified
Expiration Date	12 months from date of receipt.
Application Notes	This product is a concentrated stock solution and should be diluted appropriately with distilled, deionized water or equivalent to its final working concentration. This buffer contains untitrated 5.0 M Potassium Acetate. Meticulously prepared using ultra pure reagents dissolved in highly polished pharmaceutical grade deionized water treated with diethyl pyrocarbonate (DEPC).