

## Product Datasheet

### Tm Buffer orb348576

<b>Description</b>	Tm Buffer
<b>Hazard Information</b>	Non-Toxic
<b>Target</b>	Tm Buffer
<b>Form/Appearance</b>	Liquid (sterile filtered)
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only.
<b>Purity</b>	This product contains highly purified
<b>Expiration Date</b>	12 months from date of receipt.
<b>Application Notes</b>	This product is ready-to-use as a working 1X solution and requires no further dilution. This buffer consists of 50 mM Tris HCl, 10 mM Magnesium Sulfate at pH 7.5 and was prepared using ultra pure reagents dissolved in pharmaceutical grade deionized water.