

Product Datasheet

Hdac Assay Buffer orb348591

Description	Hdac Assay Buffer
Hazard Information	Non-Toxic
Target	Hdac Assay Buffer
Form/Appearance	Liquid (sterile filtered)
Concentration	10X
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Purity	This product contains highly purified
Expiration Date	12 months from date of receipt.
Application Notes	This product is a 10X concentrated stock solution. Prepare a 1X working solution by diluting 1 part 10X concentrate with 9 parts distilled-deionized water or equivalent. Diluted (1X) HDAC Assay buffer consists of 25 mM Tris Cl, 137 mM NaCl, 2.7 mM KCl and 1.0 mM MgCl ₂ at pH 8.0.