

Product Datasheet

Avidin Texas Red Conjugated orb348742

Description	Avidin Texas Red Conjugated
Tested Applications	ICC, IF
Hazard Information	Non-Toxic
Target	Avidin Texas Red
Form/Appearance	Lyophilized: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Purity	Avidin Texas Red Conjugated was prepared from chromatographically purified avidin from egg white followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Avidin. No rea
Expiration Date	12 months from date of receipt.
Application Notes	Avidin Texas Red Conjugated is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting.