

Product Datasheet

Beta galactosidase orb348649

Description	Onpg Elisa Beta-Galactosidase Substrate
Tested Applications	ELISA
Hazard Information	Non-Toxic
Target	Beta galactosidase
Form/Appearance	Solid
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Purity	Purity: > 95% Melting Range: 174° to 177° C Molecular Weight: 301.25 Formula: C ₁₂ H ₁₅ NO ₈ Grade: UltraPure Special Precautions: Use good judgement and appropriate handling procedures when weighing.
Expiration Date	12 months from date of receipt.
Application Notes	ONPG ELISA Beta-Galactosidase Substrate is a useful detection for ELISA.