

Product Datasheet

1X TBS pH 7.8 orb420109

Description	1X TBS pH 7.8
Hazard Information	Non-Toxic
Target	1X TBS pH 7.8
Form/Appearance	Liquid (sterile filtered)
Concentration	1X
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Expiration Date	12 months from date of receipt.
Application Notes	This product is a 1X working solution. 1X Tris Buffered Saline consists of 0.1M Tris HCl and 0.15 M Sodium Chloride, pH 7.8, meticulously prepared using ultra pure reagents dissolved in pharmaceutical grade water (WFI).