

## Product Datasheet

### Cell Lysis Buffer orb420115

<b>Description</b>	Cell Lysis Buffer
<b>Hazard Information</b>	Non-Toxic
<b>Target</b>	Cell Lysis Buffer
<b>Form/Appearance</b>	Liquid (sterile filtered)
<b>Concentration</b>	1X
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only.
<b>Expiration Date</b>	12 months from date of receipt.
<b>Application Notes</b>	The buffer is ready to use for all types of human or animal cells. Lysed cells can be used for Western Blot, co-Immunoprecipitation, Chromatin IP (ChIP), enzyme activity assays, and various other cell based assays.