

## Product Datasheet

### DNA Lambda Phage orb420125

<b>Description</b>	DNA Lambda Phage
<b>Hazard Information</b>	Non-Toxic
<b>Target</b>	DNA Lambda Phage
<b>Form/Appearance</b>	Liquid (sterile filtered)
<b>Concentration</b>	438 µg/ml
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only.
<b>Expiration Date</b>	12 months from date of receipt.
<b>Application Notes</b>	A solution in 10 mM Tris-Cl and 1 mM EDTA, pH 8.0. Note: A higher sample load may be required to clearly see the 564 and 125 base pair bands.)