

Product Datasheet

DNA Lambda Phage orb420125

Description	DNA Lambda Phage
Hazard Information	Non-Toxic
Target	DNA Lambda Phage
Form/Appearance	Liquid (sterile filtered)
Concentration	438 µg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Expiration Date	12 months from date of receipt.
Application Notes	A solution in 10 mM Tris-Cl and 1 mM EDTA, pH 8.0. Note: A higher sample load may be required to clearly see the 564 and 125 base pair bands.)