

Product Datasheet

Proteinase K (Tritirachium album limber) orb420129

Description	Proteinase K (Tritirachium album limber)
Hazard Information	Non-Toxic
Target	Proteinase K (Tritirachium album limber)
Form/Appearance	Lyophilized
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Expiration Date	12 months from date of receipt.
Application Notes	Proteinase K is specifically formulated for molecular cloning, nucleic acid research and protein/peptide structural analysis. One unit liberates one μ mole of Folin positive amino acids, measured as tyrosine, at 37° C at pH 7.5 using urea denatured hemoglobin as the substrate.