

## Product Datasheet

### Ribonuclease A orb420130

<b>Description</b>	Ribonuclease A
<b>Hazard Information</b>	Non-Toxic
<b>Target</b>	Ribonuclease A
<b>Form/Appearance</b>	Liquid (sterile filtered)
<b>Concentration</b>	5 mg/mL
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only.
<b>Expiration Date</b>	12 months from date of receipt.
<b>Application Notes</b>	One unit will cause an increase in the absorbance of 1.0 OD units at 260 nm at 37° C and pH 5.0 when yeast ribosomal RNA is hydrolyzed to acid soluble oligonucleotides.