

Product Datasheet

Protease *S. aureus* orb420132

| | |
|---------------------------|---|
| Description | Protease <i>S. aureus</i> |
| Hazard Information | Non-Toxic |
| Target | Protease <i>S. aureus</i> |
| Form/Appearance | Lyophilized |
| Concentration | 1 mg/mL |
| Storage | Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles. |
| Note | For research use only. |
| Expiration Date | 12 months from date of receipt. |
| Application Notes | Protease <i>S. aureus</i> has shown that one unit causes a change of 0.001 A280 nm per minute at 37° C at pH 7.8 using casein as a substrate. Specific conditions for reactivity should be optimized by end user. |