

Product Datasheet

Guinea Pig Serum orb436694

Description	Guinea Pig Serum
Reactivity	Guinea pig
Tested Applications	Functional Assay, IHC-Fr, IHC-P
Hazard Information	Non-Toxic
Target	Guinea Pig Serum
Form/Appearance	Liquid
Storage	Store at 4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Note	For research use only.
Purity	Serum
Expiration Date	12 months from date of receipt.