

Product Datasheet

Elisa Igm Reducing Assay Diluent orb436764

Description	Elisa Igm Reducing Assay Diluent
Tested Applications	ELISA
Hazard Information	Non-Toxic
Target	Elisa Igm Reducing Assay Diluent
Form/Appearance	Ready To Use Phosphate buffered saline
Storage	Store at 4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Note	For research use only.
Purity	Ready to use
Expiration Date	12 months from date of receipt.