

---

## Product Datasheet

### Red Cell Lysing Buffer orb436774

<b>Description</b>	Red Cell Lysing Buffer
<b>Tested Applications</b>	FACS
<b>Hazard Information</b>	Non-Toxic
<b>Target</b>	Red Cell Lysing Buffer
<b>Form/Appearance</b>	Reagent
<b>Storage</b>	Store at 4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Note</b>	For research use only.
<b>Expiration Date</b>	12 months from date of receipt.