

Product Datasheet

Antigen Retrieval Buffer, Universal orb436807

Description	Antigen Retrieval Buffer, Universal
Tested Applications	IHC-P
Hazard Information	Non-Toxic
Target	Antigen Retrieval Buffer, Universal
Form/Appearance	3 x Concentrate Biorbyt's Antigen unmasking fluid is supplied as 3x concentrated fluid. Crystals in the fluid will dissolve during heating before use.
Storage	Store at 4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Note	For research use only.
Expiration Date	12 months from date of receipt.