

Product Datasheet

Human GCSF ELISA Kit orb50023

Description	Sandwich High Sensitivity ELISA kit for Quantitative Detection of Human GCSF ELISA Kit
Reactivity	Human
Range	31.2pg/ml-2000pg/ml
Sensitivity	
Tested Applications	ELISA
Target	GCSF. No detectable cross-reactivity with any other cytokine.
Storage	Store at 4°C for 6 months and at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice).
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	Typical Data Sample Data: Concentration(pg/ml) 0 31.2 62.5 125 250 500 1000 2000 O.D. 0.011 0.135 0.279 0.464 0.928 1.496 1.956 2.308 TMB Incubation Time: 25-30min