

Product Datasheet

hsa-let-7a RT-PCR Primer Set orb220643

Description	Stem-loop RT primer and PCR primer set to detect mature miRNA hsa-let-7a expression.
Reactivity	Human
Sensitivity	10 copies of a standard mimics and as little as 25 pg of total RNA
Tested Applications	PCR
Target	hsa-let-7a
Storage	Store at 4°C for short term and at -20°C for long term storage
Note	For research use only.
Expiration Date	12 months from date of receipt.
Application Notes	The miRNAs RT-PCR primer set contains a stem-loop RT primer and a pair of specific PCR primer set. Combined with the customs RT and PCR reagent. It provides a sensitive and specific method of real-time PCR for quantificating of miRNAs from total RNA samples. It is Ideal for validating microRNA expression profiling data and high-throughput applications. Trait lies in even highly homogized miRNAs can be accurately dicriminated ensured by stemloop RT primer and high specific miRNAs primer set.