

Product Datasheet

ZSCAN1 antibody orb165171

Description	Rabbit polyclonal antibody to ZSCAN1 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide between 335-362 amino acids from the C-terminal region of human ZSCAN1.
Target	ZSCAN1
Form/Appearance	Supplied in PBS with 0.09% (W/V) sodium azide
Storage	Store at 2-8°C for short term and long term at -20°C storage, Small aliquots to prevent freeze-thaw cycles
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Peptide Affinity purified
Uniprot ID	Q8NBB4
NCBI	NP_872378.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

