

## Product Datasheet

### TXNRD3 antibody orb166341

<b>Description</b>	Rabbit polyclonal antibody to TXNRD3
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	This TXNRD3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 304-333 amino acids from the Central region of human TXNRD3.
<b>Target</b>	TXNRD3
<b>Form/Appearance</b>	PBS with 0.09% (W/V) sodium azide
<b>Storage</b>	Store at 2-8°C for up to 6 months. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Peptide Affinity purified
<b>MW</b>	70.683 kDa
<b>Uniprot ID</b>	<a href="#">Q86VQ6</a> (Human)
<b>NCBI</b>	<a href="#">NP_001166984.1</a>
<b>Entrez</b>	<a href="#">114112</a>
<b>Dilution Range</b>	WB: 1:1000
<b>Expiration Date</b>	12 months from date of receipt.

