

Product Datasheet

SQSTM1 (phospho-S207) antibody orb166801

Description	Rabbit polyclonal antibody to SQSTM1 (phospho-S207).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT, ELISA
Immunogen	This SQSTM1 Antibody is generated from rabbits immunized with a KLH conjugated synthetic phosphopeptide corresponding to amino acid residues surrounding S207 of human SQSTM1.
Target	SQSTM1 (phospho-S207)
Form/Appearance	PBS with 0.09% (W/V) sodium azide
Storage	Store at 2-8°C for up to 6 months. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Peptide Affinity purified
MW	47.687 kDa
Uniprot ID	Q13501 (Human)
NCBI	NP_001135770.1
Entrez	8878
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

