

Product Datasheet

Sqstm1 antibody orb305940

Description	Rabbit polyclonal antibody to Sqstm1
Species/Host	Rabbit
Reactivity	Human, Mouse
Predicted Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	This Phospho-Sqstm1(S351) antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 346-379 amino acids from the human region of human Phospho-Sqstm1(S351).
Target	Sqstm1. This antibody has specificity with Sqstm1
Form/Appearance	PBS with 0.09% (W/V) sodium azide
Storage	Store at 2-8°C for up to 6 months. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified
MW	48.163 kDa
Uniprot ID	Q64337 (Human)
Entrez	18412
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

