

## Product Datasheet

### Pbs Fish Gel Concentrate 10X orb348587

<b>Description</b>	Pbs Fish Gel Concentrate 10X
<b>Tested Applications</b>	ELISA, WB
<b>Hazard Information</b>	Non-Toxic
<b>Target</b>	Pbs Fish Gel Concentrate 10X
<b>Form/Appearance</b>	Liquid (sterile filtered)
<b>Concentration</b>	10X
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only.
<b>Purity</b>	This product contains highly purified
<b>Expiration Date</b>	12 months from date of receipt.
<b>Application Notes</b>	This product is a 10X concentrated stock solution. Prepare a 1X working solution by diluting 1 part 10X concentrate with 9 parts distilled-deionized water or equivalent. 10X PBS Fish Gel Concentrate consists of 0.2M potassium phosphate, 1.5M sodium chloride and Fish Gelatin at pH 7.2. A proprietary combination of stabilizers and preservatives are used that are azide and mercury free.