

Product Datasheet

14-3-3 antibody orb624373

Description	Rabbit polyclonal antibody to 14-3-3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, IP, WB
Immunogen	tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, theta polypeptide
Target	14-3-3
Form/Appearance	Liquid: PBS with 0.02% sodium azide and 50% glycerol pH 7.3
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles(Shipped with wet ice.)
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	>95% as determined by SDS-PAGE, Immunogen affinity purified
MW	28 kDa
Uniprot ID	P27348
Dilution Range	WB: 1:1000-1:10000, IP: 1:200-1:2000, IHC-P: 1:20-1:200, IF: 1:20-1:200
Expiration Date	12 months from date of receipt.